



# KN95 Mask

## User Instructions

Failure to follow the instructions and warnings on the use of this mask during all times of exposure can reduce the effectiveness of the mask and could result in illness or disability.

As with any respiratory device, the wearer must be adequately trained prior to use.

**Before use, the wearer should always check the mask is in good condition – no dirt, damage etc.**

### Fitting Instructions:

1. Hold the mask in position over the nose and mouth.
2. Pull the head harness behind the ears, attached the head harness to the retaining clip to improve comfort and prevent leakage.
3. Ensure the nose clip is securely moulded around the nose, resting the ends against the cheek to obtain a good seal.
4. To check for a proper fit, cup both hands over the mask and exhale vigorously. If air leaks around the nose, tighten the nose clip, if air leaks around the edge, reposition the head harness for a better fit.
5. Repeat adjustments until the mask is sealed properly.
6. If a proper seal cannot be achieved do not use.

### Warnings:

1. This mask shall not be used for more than one shift
2. Never substitute, modify, add or omit parts in the configuration as specified by the manufacture.
3. This mask helps protect against certain particulate contaminants but does not completely eliminate exposure to the risk of the contracting disease or infection.
4. It is unlikely that the requirements for leakage will be achieved if facial hair passes under the face seal.
5. Change the mask immediately if breathing becomes difficult or the mask becomes damaged or distorted.
6. Change the mask if a proper face seal cannot be achieved.

Do not use the respirator or enter or stay in a contaminated area under the following circumstances:

- Atmosphere contains less than 19.5% oxygen
- If you smell or taste contaminant
- For protection against gasses or vapors
- Contaminants or their concentrations are unknown or immediately dangerous to life or health
- For sandblasting, paint-spraying operations and asbestos treatment
- In explosive atmospheres