

# bastion


DISPOSABLE GLOVES

Progenics™ Glove Dispenser System



bastion

Progenics™ Cuff First



- Reduces risk of contamination.
- Gloves dispensed individually.
- Cost effective due to reduced waste.
- Gloves kept away from contaminated surfaces.

RE-ORDER CODE:

The Progenics Cuff-First dispensing system helps to reduce glove wastage and reduce costs. As only one glove is dispensed at a time, wastage of gloves and costs is reduced considerably. Studies in Norway and the USA have shown a reduction in glove consumption of between 20% and 40% when dispensed using a cuff-first dispensing system.

bastion

Ultimate solutions in barrier protection

T 0800 864 725 F 0800 864 725 E [sales@bastionpacific.co.nz](mailto:sales@bastionpacific.co.nz)  
W [www.bastionpacific.co.nz](http://www.bastionpacific.co.nz)

# Progenics™ Glove Dispenser System

Recent studies show gloves dispensed from Bastion Progenics™ System have **92.5% less bacterial contamination** than traditionally dispensed gloves.



## Progenics™ System

- Gloves are dispensed cuff first – reducing the risk of contamination by not touching the finger and palm areas of the glove.
- Fast and easy to use.
- Cost effective due to reduced waste as gloves are dispensed individually.
- Unused gloves cannot be reinserted into the box further reducing cross contamination.
- The acrylic wall mounted dispenser bracket keeps gloves away from contaminated surfaces and boxes are easily removed and replaced.



## Progenics Dispenser Brackets

- Acrylic wall mounted dispenser bracket keeps gloves away from contaminated surfaces.
- Boxes are easily removed and replaced.
- Brackets are supplied in kits with screw fixings, double sided tape and mounting instructions.



Item Code	Size
1000/113	Single Wall
1000/117	Triple Wall



## Progenics Vinyl P/F Clear Gloves

Item Code	Size
100/611	Small (8cm)
100/612	Medium (9cm)
100/613	Large (10cm)
100/614	X-Large (11cm)

## Progenics Vinyl P/F Blue Gloves

Item Code	Size
100/621	Small (8cm)
100/622	Medium (9cm)
100/623	Large (10cm)
100/624	X-Large (11cm)



## Progenics Nitrile P/F Blue Gloves

Item Code	Size
100/631	Small (8cm)
100/632	Medium (9cm)
100/633	Large (10cm)
100/634	X-Large (11cm)



Patent No.: AU 2009284021

# bastion

Ultimate solutions in barrier protection

T 0800 864 725 F 0800 864 725 E sales@bastionpacific.co.nz  
W www.bastionpacific.co.nz